

1m4 AURISclear

Professional ear irrigation system for ear cleaning

- ▶ For gentle cleaning of the ear canal with water
- ▶ With integrated thermometer for temperature control
- ▶ AURISlips 3-stream disposable tips with Luer-Lock connection
- ▶ Content: 1 x 500 ml spray bottle, 1 x tip adapter, 1 x drip tray and 20 x AURISlips disposable tips



Technical specifications

Material:	Spray bottle: Polyethylene (HDPE/PE), Tip: Thermoplastic Elastomer (TPE) Drip tray and adapter: Acrylonitrile Butadiene Styrene (ABS)
Dimension:	280 x 140 x 90 mm
Weight:	130 g
HS-Code*:	90189084 DE

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
AURISclear	1 pieces	10 x 1 piece	4745010564321	10029.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland