

1m4 AURISsingle

Ear irrigation syringe for safe removal of earwax

- ▶ For gentle cleaning of the ear canal with water
- ▶ Safe application without risk to the eardrum
- ▶ Contents: 1 x 20 ml syringe, 1 x three-stream attachment (AURISstip)



Technical specifications

Material:	Polypropylene (PP), Thermoplastic Elastomer (TPE)
Dimension:	50 x 40 x 130 mm
Weight:	30 g
HS-Code*:	90183110 DE

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
AURISsingle	1 piece	100 x 1 piece	4745010564499	10040.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland