

1m4 AURISlips

Disposable tips for ear irrigation systems for ear cleaning

- ▶ 3-stream disposable tips with Luer-Lock connection
- ▶ For gentle cleaning of the ear canal with water
- ▶ for use with: AURISclear, devices using OtoClear®** by Bionx®** and Luer-Lock syringes



Technical specifications

Material:	Thermoplastic Elastomer (TPE)
Dimension:	160 x 110 x 50 mm
Weight:	100 g
HS-Code*:	90189084 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
AURISlips	40 pieces	32 x 40 pieces	4745010564321	10029.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland

** OtoClear® and Bionx® are registered trademarks of BIONIX® LLC, Maumee, Ohio, USA, and are not affiliated with 1M4 OÜ in any way.