


1m4 ECG-Ruler

For evaluation and interpretation of ECGs

- ▶ Handy interpretation aid for precise EKG evaluation
- ▶ PQ-, QRS-, QT-, scale for 25 mm/s + 50 mm/s and amplitudes
- ▶ Detailed, clear representation
- ▶ Easy to clean and abrasion-resistant
- ▶ Custom branding available



Technical specifications

Material: Polyvinyl chloride (PVC)

Dimension: 242 x 50 mm

Weight: 16 g

HS-Code*: 90172090 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-N0.
EKG-Lineal	1 piece	100 x 1 piece	4745010564314	10027.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland