

# 1m4 ECG-Caliper

For evaluation and measurement of ECGs

- ▶ For manual measurement of ECG amplitudes
- ▶ No unintended change of the set measurement distance
- ▶ Resistance can be adjusted via screws
- ▶ Cover with ruler
- ▶ Easy to clean
- ▶ Custom branding available



## Technical specifications

Material:	Brass (chrome-plated), Acrylonitrile Butadiene Styrene Copolymer (ABS)
Dimension:	115 mm
Weight:	38 g
HS-Code*:	90178090 CN

## Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
ECG-Caliper	1 piece	25 x 1 piece	4745010564307	10028.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: [info@medichem-online.de](mailto:info@medichem-online.de)



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland