

1m4 Instrument Brush

For cleaning instruments in three variants

- ▶ For effective cleaning of instruments
- ▶ Also removes persistent dirt
- ▶ Variants with bristles made of: brass, nylon, or stainless steel
- ▶ With ergonomic handle
- ▶ Autoclavable, heat resistant up to 134 °C



Caution: When using cleaning brushes with metal bristles, please handle with care, as the bristles can damage sensitive surfaces.

Technical specifications

Material: Polypropylene (PP), Brass, Nylon, Stainless Steel

Dimension: 180 mm, bristle length: 13 mm

Weight: *22 g, **18 g

HS-Code*: 96039099 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Brass*	1 piece	10 x 1 piece	4745010564451	12002.01
Nylon**	1 piece	10 x 1 piece	4745010564468	12002.02
Stainless Steel*	1 piece	10 x 1 piece	4745010564475	12002.03

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de