


1m4 LH-Test

Vision Test for Diagnosing Visual Acuity

- ▶ Specifically designed for children and individuals with communication challenges
- ▶ Standardized and versatile – ideal for routine check-ups such as developmental assessments, as well as for school, occupational, and neurological evaluations.
- ▶ Results can be documented in various formats: Decimal Visual Acuity (-0.3 to 1.0), LogMAR (0.1 to 2.0), M-Value (1.5 to 30), and Snellen Fraction (3/6 m)
- ▶ Practical ring binder containing 22 symbol charts, occluder card, and response card
- ▶ Evaluation sheet available in PDF format



Technical specifications

Material:	Paper, Steel, Polyethylene (PE)
Dimension:	150 x 148 x 5 mm
Weight:	132 g
HS-Code*:	49019900 DE

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
LH-Test	1 piece	10 x 1 piece	4745010564598	10045.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland