

1m4 Tape Measure

For easy measurement of circumferences

- ▶ Robust housing with retractable mechanism
- ▶ Continuously extendable tape measure
- ▶ Scale in cm and inches
- ▶ Measurement range: 0-203 cm, graduation 1 mm;
0-80 inches, graduation 1/16 inch
- ▶ Custom branding available



Technical specifications

Material:	Acrylonitrile Butadiene Styrene (ABS), Polyvinyl Chloride (PVC)
Dimension:	Ø 70 x 20 mm
Weight:	50 g
HS-Code*:	90178010 DE

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Tape Measure	1 piece	25 x 1 piece	4745010564543	10041.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany
Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de