


1m4 Nose Clip Pads

Pads for nasal clips for spirometric examination

- ▶ For use with nose clips from various manufacturers
- ▶ Soft padding
- ▶ For single use



Technical specifications

Material:	Ethylene-Vinyl-Acetate-Copolymer (EVAC)
Dimension:	10 x Ø 17 mm
Weight:	0,5 g
HS-Code*:	392690979 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Nose Clip Pads	100 pieces	50 x 100 pieces	4745010564222	10020.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland