

1m4 Scoliometer

Diagnostic instrument for measuring the spine

- ▶ For measuring the lateral curvature of the spine (trunk asymmetry)
- ▶ With notch for spinous processes
- ▶ Applications: Pediatrics, Orthopedics, and Physiotherapy
- ▶ Easy to clean
- ▶ Plastic free packaging
- ▶ Custom branding available



MD

Technical specifications

Material:	Acrylonitrile-Butadiene-Styrene-Copolymer (ABS), Polymethyl methacrylate (PMMA)
Dimension:	168 x 54 x 10 mm
Weight:	60 g
HS-Code*:	90318020 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Scoliometer	1 piece	25 x 1 piece	4745010564291	10026.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de

CH REP

TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland