


1m4 Goniometer

For measuring joint angles

- ▶ For measuring joint position and mobility according to the neutral-zero method
- ▶ Suitable for orthopedics, physiotherapy, and rehab medicine
- ▶ With printed degree scale and ruler (graduation in centimeters and inches)
- ▶ Robust, flexible, and transparent plastic
- ▶ Easy to clean
- ▶ Custom branding available



MD

Technical specifications

Material: Polyvinyl chloride (PVC)

Dimension: 215 x 45 mm

Weight: 25 g

HS-Code*: 90172090 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Goniometer	1 piece	50 x 1 piece	4745010564147	10014.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de

CH REP

TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland